




Across-Visit Co-Registration of T1 Images





Across-Visit Co-Registration of T1 Images

- Visit 1 files are in the visit 1 T1 space and visit 2 files are in the visit 2 T1 space
- \Box We will register the follow-up/visit 2 T1 scan to the visit 1/baseline scan
- □ Next slide:
 - □ take the skull-stripped T1 from visit 2, register it to the skull-stripped T1 from visit 1
 - use this transformation to put the skull-stripped T2 and skull-stripped FLAIR from visit 2 into the T1 space from visit 1





Across-Visit Co-Registration Skull-Off Images



Plotting Registered T1 Images Side-By-Side


```
ss_t1=masked_imgs[[1]]
visit_2_t1=readNIfTI(outfiles2[1],reorient=FALSE)
double_ortho(ss_t1,visit_2_t1)
```







Plotting Registered T1 Image Overlay

ortho2(ss_t1, visit_2_t1, col.y=alpha(hotmetal(), 0.25))






Across-Visit Co-Registration Skull-On Images



T1 Skull-On Images Side-By-Side


```
t1_skull=readNIfTI(outfiles[1],reorient=FALSE)
v2_t1_skull=readNIfTI(outfiles2_skull[1],reorient=FALSE)
double_ortho(t1_skull, v2_t1_skull)
```






T1 Skull-On Images, Overlay

ortho2(t1_skull, v2_t1_skull,col.y=alpha(hotmetal(),0.25))





Overview

- □ Registration within a subject can be done in R
 - ants_regwrite wraps around the reading/writing of images and applying transformations
 - double_ortho and ortho2 can provide some basic visual checks to assess registration quality
- Once images are registered in the same space,
 operations can be applied to all the images:
 - Masking with a brain mask
 - Transforming images to new spaces with one modality